

BIOCLEAR

IVD In vitro diagnostic medical device



Xylene for histology INSTRUCTIONS FOR USE

REF Catalog number: BC-1L (1000 ml)

BC-5L (5000 mL)

BC-10L (10000 mL)

Introduction

Xylene is the most commonly used agent for clearing the microscopic samples. It is done as a part of the sample processing procedure and slide staining procedure. Clearing agents are used as an intermedium for mixing alcohol with embedding paraffin (BioWax) and resin alcohol/covering gel. They should be soluble in dehydrating and infiltrating medium. BioGnost's BioClear enables high quality and effective sample/slide clearing and it provides superior results. Adequately dehydrated and cleared tissue sample becomes transparent and appropriate for microscopic analysis of the slide.

Product description

- **BIOCLEAR** - Xylene for clearing, easily miscible with alcohol solutions, paraffin and synthetic covering/mounting media.

Product use

- BioClear is a clearing agent and it is used as an intermedium for mixing alcohol with the paraffin embedding wax and alcohol with resins/covering gels.
- BioClear is used in sample processing procedure and staining slides, microscope sections, cytologic and histopathologic smears.
- As a part of sample processing procedure, BioClear is used to remove alcohol from the sample and infiltrate/immerse the sample with paraffin (BioWax PLUS 56/58).

Note: the sample is treated with BioClear mixed with alcohol and it is gradually being removed from the sample, while BioClear remains in it. The sample is then treated with paraffin. It gradually removes BioClear in a similar way and infiltrates/immerses the sample.

- As a part of the slide staining procedure, BioClear is used for removing the paraffin from the slide and infiltration/immersion of the slide using alcohol (Histanol 70, Histanol 80, Histanol 95, Histanol 100).

Recommended procedure for removing paraffin from the section (deparaffination):

Deparaffine in xylene (BioClear)	3 exchanges, 2 min each
----------------------------------	-------------------------

Note: the sample is treated with BioClear mixed with paraffin and it is gradually being removed from the sample, while BioClear New remains in it. The slide is then treated using alcohol (series of ascending concentrations) that removes BioClear from the slide in the same way and dehydrates the slide.

- After staining clear the sections using BioClear medium; the section must be covered afterwards.

Recommended procedure for section clearing:

Clear in xylene (BioClear)	2 exchanges, 2 min each
----------------------------	-------------------------

- Cover/mount the samples treated with BioClear using xylene-based media, such as BioMount, BioMount High, BioMount M, BioMount DPX, BioMount DPX High, BioMount DPX Low, BioMount DPX New, or BioMount C or an aliphatic hydrocarbon based mounting medium, such as BioMount New and BioMount New Low.

Result

The result of sample clearing (part of the sample processing procedure) is a simplified, high quality and complete infiltration/immersing the sample with paraffin, enabling easy and high quality cutting the sample using microtome and cutting it into strips.

The result of slide clearing (part of the sample staining) is a simplified and high quality sample dehydration using series of alcohols of ascending concentrations, enabling high quality sample covering and preparing the slide for analysis.

Note

Depending on personal preferences and standard operating laboratory procedures, period of clearing using BioClear may vary. BioClear may cause sample/slide damage if the exposure lasts for too long.

Preparing the sample and diagnostics

Use only appropriate instruments for collecting and preparing the samples. All the samples must be processed with the most modern technology and be visibly marked. Follow the manufacturer's instructions for handling. In order to avoid mistakes, staining must be conducted by a trained professional. Only trained medical personnel may make a diagnosis. Use only microscope according to standards of the medical diagnostic laboratory. In order to avoid an erroneous result, a positive and negative check is advised before application.

Safety at work and environmental protection

Handle the product in accordance with safety at work and environmental protection guidelines. Used solutions and out of date solutions should be disposed of as special waste in accordance with national guidelines. Chemicals used in this procedure could pose danger to human health. Tested tissue specimens are potentially infectious. Necessary safety measures for protecting human health should be taken in accordance with good laboratory practice. Act in accordance with signs and warnings notices printed on the product's label, as well as in BioGnost's material safety data sheet which is available on demand.

Storing, stability and expiry date

Keep BioClear in a tightly closed original package at temperature between +15°C and +25°C. Do not freeze and avoid exposing to direct sunlight. Date of manufacture and expiry date are printed on the product's label.

References

1. Carson, F. L. (2007), Histotechnology, 2nd ed. Singapore
2. Cook, D. J. (2006): Cellular pathology, 2nd ed. Banbury: Scion Publishing Ltd.
3. Kiernan, J. A. (2008) Histological and histochemical methods, 4th ed. Bloxham: Scion Publishing Ltd.

BC-X, V13-EN7, 02 July 2019, AK/IŠP

Refer to the supplied documentation	Storage temperature range	Number of tests in package	Product code	European Conformity
Refer to supplied instructions	Keep away from heat and sunlight	Valid until	Lot number	Manufacturer
For <i>in vitro</i> diagnostic use only	Keep in dry place	Caution - fragile		

BIOGNOST Ltd.
Medjugorska 59
10040 Zagreb
CROATIA
www.biognost.com

