# **BIOGNOST**®

# **HEMATOXYLIN INSTANT KIT**

CE

IVD In vitro diagnostic medical device

# For fast and simple preparation of 1 liter of Hematoxylin Instant, used in histology and cytology

# **Contains two components**

## **INSTRUCTIONS FOR USE**

REF Catalogue number: HEMI-K-1L (for preparing 1 L of Hematoxylin Instant) HEMI-K-6L (for preparing 6 x 1 L of Hematoxylin Instant)

#### Introduction:

BioGnost's Hematoxylin Instant kit is specially prepared for fast and practical preparation of Hematoxylin Instant that can be used for staining histology and cytology samples. It contains two specially stabilized components that (by dissolving in distilled water) create a formulation of hematoxylin that may be used as a substitute for Harris' or Gill's hematoxylin. Hematoxylin Instant provides very precise and clear results of cellular nuclei, and may also be used in progressive and in regressive staining methods. It also does not contain oxidant based on mercury chloride. That makes it environment-friendly.

#### **Product description:**

HEMATOXYLIN INSTANT KIT - kit for preparing Hematoxylin Instant for nuclear staining in histopathology and cytology.

Hematoxylin Instant kit contains:	for preparing 1 L of Hematoxylin Instant (HEMI-K-1L)	for preparing 6 x 1 L of Hematoxylin Instant (HEMI-K-6L)				
Hematoxylin Instant - component A	1 x HEMI-KA	6 x HEMI-KA				
Hematoxylin Instant - component B	1 x HEMI-KB	6 x HEMI-KB				

#### Preparing Hematoxylin Instant

1.	Pour all the contents of components A and B of Hematoxylin Instant in the container that can fit more than 1 L of fluid.								
2	Add 1 L of distilled (demi) water. DO NOT USE TAP WATER. If Hematoxylin Instant is used for staining cytology sections or in progressive								
۷.	staining methods, dissolve components A and B in 1.5 do 2 L of distilled (demi) water.								
	Note: in order to achieve better stability, less precipitate and longer expiry date, we recommend 900 mL of distilled water and 100 mL of ethylene glycol								
	(CAS: 107-21-1). For cytology or progressive staining - dissolve the components in 1350 mL of distilled water and 150 mL of ethylene glycol or, for								
	even more specific staining, in 1800 mL of distilled water and 200 mL of ethylene glycol.								
3.	Mix the solution well until the contents melt. Wait for at least 60 minutes before first staining.								
	Note: for optimal results wait 12 hours before first staining.								
4.	Filter Hematoxylin Instant before use. Precipitate and metallic surface are normal occurrence for this type of product.								
5	Keep Hematoxylin Instant in a tightly closed original package at temperature between +15°C and +25°C. Keep in dry places, do not freeze and								
J.	avoid exposing to direct sunlight. Shelf life of prepared Hematoxylin Instant is 6-9 months.								

#### Using Hematoxylin Instant:

#### a) in histological sections staining procedure

#### Preparing histological sections for staining

- Fix the tissue sample tightly (4% NB Formaldehyde, 10% NB Formaldehyde), rinse with water and dehydrate through series of ascending alcohol solutions (Histanol 70, Histanol 80, Histanol 95 and Histanol 100).
- Clear the sample with intermedium; in xylene (BioClear) or in a xylene substitute (BioClear New).
- Infiltrate and fit the sample in paraffin (BioWax Plus 56/58, BioWax 52/54, BioWax 56/58, BioWax Blue, BioWax Micro).
- Cut the paraffin block to 4-6  $\mu$ m slices and place them on a VitroGnost glass slide.

#### Hematoxylin and eosin (HE) staining procedure, regressive

Deparaffinize the section in xylene (BioClear) or in a xylene substitute (BioClear New)	3 exchanges, 2 min each				
Rehydrate in descending series of alcohols (Histanol 100, Histanol 95) and in distilled (demi) water	10 dips in each of the 3 exchanges				
Stain using Hematoxylin Instant	2-5 minutes				
Note: incubation period is shorter with progressive staining	30 seconds to 2 min				
Rinse with distilled (demi) water	1 min				
Differentiate in Acid alcohol	3-10 fast dips				
Note: This step removes excessive hematoxylin. Discoloration of the nuclei can occur if the section is treated with the differentiation agent for too long. Skip this step with progressive staining					
Rinse in distilled (demi) water					
Blue using Scott's solution or Bluing reagent	1 min				
Note: If the mentioned reagents are not available, the section should be blued using indirect stream of water	3-5 minutes				
Use Eosin Instant, alcoholic or aquoues until the section is optimally stained. Eosin Instant, alcoholic will achieve optimal staining for 30 seconds - 1 minute and Eosin Instant, aqueous in under 3	30 seconds - 3 min				
minutes.					
Note: Eosin Instant may be substituted with any BioGnost's aqueous or alcoholic Eosin solution					
Rinse under tap water	30 seconds				
Dehydrate using 95% alcohol (Histanol 95)	10-15 dips in each of the 2				
	exchanges				
Completely rehydrate by using 100% alcohol (Histanol 100)	10-15 dips in each of the 3				
	exchanges				
Clear the section in xylene (BioClear) or in a xylene substitute (BioClear New)	2 minutes in each of the 2 exchanges				
	Deparaffinize the section in xylene (BioClear) or in a xylene substitute (BioClear New)   Rehydrate in descending series of alcohols (Histanol 100, Histanol 95) and in distilled (demi) water   Stain using Hematoxylin Instant   Note: incubation period is shorter with progressive staining   Rinse with distilled (demi) water   Differentiate in Acid alcohol   Note: This step removes excessive hematoxylin. Discoloration of the nuclei can occur if the section is treated with the differentiation agent for too long. Skip this step with progressive staining   Rinse in distilled (demi) water   Blue using Scott's solution or Bluing reagent   Note: If the mentioned reagents are not available, the section is optimally stained. Eosin Instant, alcoholic will achieve optimal staining for 30 seconds - 1 minute, and Eosin Instant, aqueous in under 3 minutes.   Note: Eosin Instant may be substituted with any BioGnost's aqueous or alcoholic Eosin solution   Rinse under tap water   Dehydrate using 95% alcohol (Histanol 95)   Completely rehydrate by using 100% alcohol (Histanol 100)   Clear the section in xylene (BioClear) or in a xylene substitute (BioClear New)				

Immediately after clearing apply an appropriate BioMount medium for covering/mounting on the section. If BioClear xylene was used, use one of BioGnost's mounting xylene-based media (BioMount, BioMount High, BioMount M, BioMount DPX, BioMount C, or universal BioMount New). If BioClear New xylene substitute was used, the appropriate covering agent is BioMount New. Cover the section with a VitroGnost cover glass.

#### **Result:**

Nuclei - blue

Cytoplasm, collagen, muscle fibers, erythrocytes - hues of pink to red

### b) in cytology (and gynecology samples) staining methods acc. to Papanicolaou, <u>regressive</u>

The first stage of staining procedure depends on the method the cytological sample was collected and fixated on the microscope slide.

If the sample is dry and previously fixated by using solution containing polyglycols (CitoSpray), it is necessary to keep it in a 95% alcohol solution (Histanol 95) for 10 minutes in order to remove polyglycols. If the section was fixated with a 95% alcohol solution (Histanol 95), ignore this step. During staining cytology samples (prepared by using the liquid based cytology method (LBC)) that contain low concentration of alcohol, rehydration by descending series of alcohol solutions is not necessary. The procedure begins by rinsing the section using distilled (demi) water and is then stained using Hematoxylin Instant.

1	Rehydrate in descending series of alcohols (Histanol 95, Histanol 80 and Histanol 70) and in distilled	6-8 dips in each of the 4				
1.	or demineralized water	exchanges				
2.	Stain using Hematoxylin Instant	4 min				
	Note: incubation period is shorter with progressive staining	30 sec				
3.	Rinse in distilled/demineralized water	6-8 dips				
4.	Differentiation using HCL Pap reagent or in 0.1% HCl solution	5-10 seconds				
	Note: This step removes excessive hematoxylin from the nucleus and cytoplasm. Discoloration of the nuclei can occur if the section is treated with the differentiation agent for too long.					
5.	Rinse in distilled water	6-8 dips				
6.	Blue using Scott's solution or Bluing reagent	1 min				
	Note: If the mentioned reagents are not available, the section should be blued using indirect stream of water	3-5 minutes				
7.	Debudrate in accending series of alcohols (Histanol 70, Histanol 80 and Histanol 95)	6-8 dips in each of the 3				
		exchanges				
8.	Stain using OG-6, Pap 2A reagent	2 min				
٥	Rinse using 05% alcohol in two exchanges (Histanol 05)	6-8 dips in each of the 2				
9.		exchanges				
10.	Stain using EA 31, Pap 3A reagent or EA 50, Pap 3B reagent	4 min				
11.	Rinse using 95% alcohol (Histanol 95)	6-8 dips				
12.	Dehydrate using 100% alcohol (Histanol 100)	6-8 dips				
13.	Dehydrate using 100% alcohol (Histanol 100)	3-5 minutes				
14.	Clear the section in xylene (BioClear) or in a xylene substitute (BioClear New)	6-8 dips				
15.	Clear the section in xylene (BioClear) or in a xylene substitute (BioClear New)	3-5 minutes				

Immediately after clearing apply an appropriate BioMount medium for covering/mounting on the section. If BioClear xylene was used, use one of BioGnost's mounting xylene-based media (BioMount, BioMount High, BioMount M, BioMount DPX, BioMount C, or universal BioMount New). If BioClear New xylene substitute was used, the appropriate covering agent is BioMount New. Cover the section with a VitroGnost cover glass.

#### Results

Nuclei - blue

Keratinized cells - yellow-orange

Superficial squamous epithelial cell, erythrocytes, nucleoli, cilia - pink-red

Cytoplasm of all the other cell types (parabasal and intermediate squamous cells, columnar cells, polymorphonuclear leukocytes, lymphocytes, histiocytes, adenocarcinomas, undifferentiated carcinoma cells) - green

#### Note

Time periods of staining processes are not entirely standardized and they approximately correspond to clinical and laboratory practical experience. Intensity of staining depends on the period of immersion in the dye. Real staining protocol depends on personal requests and priorities.

#### Preparing the sample and diagnostics

Use only appropriate instruments for collecting and preparing the samples. Process the samples with modern technology and mark them clearly. Follow the manufacturer's instructions for handling. In order to avoid mistakes, the staining procedure and diagnostics should only be conducted by authorized and qualified personnel. Use only microscope according to standards of the medical diagnostic laboratory.

### Safety at work and environmental protection

Handle the product in accordance with safety at work and environmental protection guidelines. Used solutions and out of date solutions should be disposed of as special waste in accordance with national guidelines. Chemicals used in this procedure could pose danger to human health. Tested tissue specimens are potentially infectious. Necessary safety measures for protecting human health should be taken in accordance with good laboratory practice. Act in accordance with signs and warnings notices printed on the product's label, as well as in BioGnost's material safety data sheet.

#### Storing, stability and expiry date

Keep Hematoxylin Instant kit in a tightly closed original package at temperature between  $+15^{\circ}$ C and  $+25^{\circ}$ C. Keep in dry places, do not freeze and avoid exposing to direct sunlight. Date of manufacture and expiry date are printed on the product's label.

#### References

- 1. Baker, J.R. (1962): Experiments on the action of mordants. 2. Aluminium-hematein. Q.J.Microsc. Sci. p103 493-517.
- 2. Conn, J. (1977): *Biological Stains*, 9<sup>th</sup> ed., Baltimore: Williams and Wilkens Co.
- 3. Harris, H.F. (1898): A new method of "ripening" haematoxylin. Microsc. Bull. (Philadelphia) Dec. 47.
- 4. Harris, H.F. (1900): On the rapid conversion of haematoxylin into haematein in staining reactions. J. Appl. Microsc. p3 777-780.

#### HEMI-K-X, V3-EN3, 10 July 2019, AK/IŠP

	, VO EIVO, TO OULY 20	10, / 101								_			
Â	Refer to the supplied documentation	°c-	Storage temperature range	$\Sigma$	Number of tests in package	REF	Product code	CE	European Conformity		BIOGNOST Ltd. Medjugorska 59 10040 Zagreb	C	E
[]i]	Refer to supplied instructions	歉	Keep away from heat and sunlight	8	Valid until	LOT	Lot number	***	Manufacturer		CROATIA www.biognost.com		
IVD	For in vitro diagnostic use only	Ť	Keep in dry place	ų	Caution - fragile					-			