

EOSIN B powder dye, C.I. 45400

IVD *In vitro* diagnostic medical device **CE**

Eosin Bluish, Imperial Red, Acid Red 91, BSC certified dye
For staining according to Romanowsky

INSTRUCTIONS FOR USE

REF Catalogue number: EB-P-25 (25 g)

Introduction

Histology, cytology and other related scientific disciplines study the microscopic anatomy of tissues and cells. In order to achieve a good tissue and cellular structure, the samples need to be stained in a correct manner. Eosin B powder dye can be used instead of Eosin Y WS dye in different staining methods. These two dyes are mutually completely interchangeable, but Eosin B provides barely visible bluish hue instead of yellowish provided by Eosin Y. It is often used in Romanowsky staining methods. It is also used as counterstain to hematoxylin during cytoplasmic staining.

Product description

- **EOSIN B** - Biological Stain Commission (BSC) certified powder dye for preparing the solution for staining according to Romanowsky

Other preparations and reagents used in preparing the dye solution:

- Glacial acetic acid (CH₃COOH)

Preparing the solutions for staining

0.1% aqueous eosin solution:

- Dissolve 1 g of Eosin B stain in 1 liter of distilled/demineralized water while stirring.
- Add 1.6 ml of glacial acetic acid.
- Filter before use.

Result

Nuclei - dark blue
 Cytoplasm, elastin, collagen - red
 Erythrocytes - red-orange

Note

The mentioned formulation is only one of the ways of preparing the dye solution. Depending on personal requests and standard laboratory operating procedures, the dye solution can be prepared according to other protocols.

Preparing the sample and diagnostics

Use only appropriate instruments for collecting and preparing the samples. Process the samples with modern technology and mark them clearly. Follow the manufacturer's instructions for handling. In order to avoid mistakes, the staining procedure and diagnostics should only be conducted by authorized and qualified personnel. Use only microscope according to standards of the medical diagnostic laboratory. Both positive and negative controls are recommended before applying.

Safety at work and environmental protection

Handle the product in accordance with safety at work and environmental protection guidelines. Used solutions and out of date solutions should be disposed of as special waste in accordance with national guidelines. Chemicals used in this procedure could pose danger to human health. Tested tissue specimens are potentially infectious. Necessary safety measures for protecting human health should be taken in accordance with good laboratory practice. Act in accordance with signs and warnings notices printed on the product's label, as well as in BioGnost's material safety data sheet which is available on demand.










Storing, stability and expiry date

Keep Eosin B powder dye in a tightly closed original package at temperature between 15°C and 25°C. Keep in dry places, do not freeze and avoid exposure to direct sunlight. Expiry date is stated on the product's label.

References

1. Conn, J. (1977): Biological Stains, 9th ed., Baltimore: Williams and Wilkins Co.

EOB-P-25, V2-EN1, 30 September 2015, VR/IŠP

	Refer to the supplied documentation		Storage temperature range		Number of tests in package	REF	Product code	CE	European Conformity
	Refer to supplied instructions		Keep away from heat and sunlight		Valid until	LOT	Lot number		Manufacturer
IVD	For <i>in vitro</i> diagnostic use only		Keep in dry place		Caution - fragile				

 BIOGNOST Ltd.
 Medjugorska 59
 10040 Zagreb
 CROATIA
 www.biognost.com

