

LEISHMAN'S SOLUTION

IVD In vitro diagnostic medical device



Polychromatic solution of eosin, Methylene Blue and azure dyes For staining in hematology and clinical cytology

INSTRUCTIONS FOR USE

REF Catalog number: LE-OT-100 (100 mL)

LE-OT-500 (500 mL)

LE-OT-1L (1000 mL)

Introduction

Romanowsky polychromatic dyes are standardly used in hematology and cytology. Various sorts of Romanowsky dyes (Giemsa, May-Gruenwald, Leishman, Wright, Jenner) contain different ratios of methylene bluing reagent used as the cation component (and the reagent-related thiazine dyes, such as azure B) and eosin Y as the anion component. Known Romanowsky effect is created by mutual interaction between cationic and anionic components, and it is being demonstrated by purple color. Staining intensity depends on the azure content, as well as azure to eosin Y ratio, while a few other factors affect the result of staining: working solution pH value and buffer solution, fixation method and dye exposure time.

Product description

LEISHMAN'S SOLUTION - Solution of eosin, Methylene Blue and azure in methanol.

Other sections and reagents that may be used in staining:

- Dehydrating/rehydrating agent, such as BioGnost's alcohol solutions: Histanol 70, Histanol 80, Histanol 95 and Histanol 100
- High-quality glass slides for use in histopathology and cytology, such as VitroGnost SUPER GRADE, VitroGnost COLOR or one of more than 30 models of BioGnost's VitroGnost glass slides
- Fixatives, such as BioGnost's Histanol M
- BioGnost's Immersion oil
- BioGnost's Buffer tablets, pH 6.8 or 7.2

Preparation of solutions

Filter undiluted Leishman's solution before use.

Buffer solutions

Buffer solution pH 6.8 (for hematology and cytology smears).

Buffer solution pH 7.2 (for staining hematology smears expected to contain blood parasites).

- Dissolve 1 buffer tablet in 1 liter of distilled water while stirring. Filter after dissolving.
Note: During the staining process it is possible to use pH 6.8 or pH 7.2 buffer solution or a combination of pH 6.8 and 7.2 buffer solutions. The process's results can differentiate in shift toward red or blue on the color spectrum.

Leishman's working solution for vertical staining

- Combine 30mL of Leishman's solution with 150 ml of distilled or demineralized water and with 20 ml of pH 6.8 or pH 7.2 buffer solution. Let it set for 10 min. Filter before use.

Working Leishman's solution for staining in automatic stainer

- Combine 50mL of Leishman's solution with 220 ml of distilled or demineralized water and with 30 ml of pH 6.8 or pH 7.2 buffer solution. Let it set for 10 min. Filter before use.

A1) Procedure of horizontal smear staining (on a staining rack)

1.	Dry the preparation	
2.	Place the section in the horizontal position and cover it with 1 ml of undiluted Leishman's solution	1 min
3.	Add 2 mL of Buffer solution pH 6.8 or pH 7.2 , gently stir and let it react	5 min
4.	Rinse with Buffer solution , pH 6.8 or pH 7.2 through two exchanges	2 exchanges, 1 minute each
5.	Dry the preparation	

A2) Procedure of vertical smear staining (in a staining jar)

1.	Dry the preparation	
2.	Immerse the section into non-diluted Leishman's solution	3 min
3.	Immerse the section into Leishman's working solution for vertical staining	6 min
4.	Rinse with Buffer solution , pH 6.8 or pH 7.2 through two exchanges	2 exchanges, 1 minute each
5.	Dry the slide	

A3) Procedure of staining smear in automatic stainer

1.	Dry the preparation	
2.	Immerse the section into non-diluted Leishman's solution	3 min
3.	Immerse the section into Leishman's working solution for staining in automatic stainer	6 min
4.	Rinse with Buffer solution , pH 6.8 or pH 7.2	1 min
5.	Rinse in tap water	2 min
6.	Dry the preparation	

	Result (pH 6.8)	Result (pH 7.2)
Nucleus	red-purple	purple
Lymphocyte cytoplasm	blue	light blue
Monocyte cytoplasm	grey-blue	blue-grey
Neutrophil granules	bright purple	light red-purple to pink
Eosinophil granules	dark red	red-pink to brown-red
Basophil granules	dark purple	dark blue-purple to black
Thrombocyte granules	purple	reddish-purple
Erythrocytes	reddish	reddish-grey
Blood parasites	/	red-purple nuclei

Note

Time periods of staining processes are not entirely standardized in clinical and laboratory practical experience. Time periods specified in the instruction approximately correspond to a longtime work practice with optimal results. Intensity of staining depends on the period of immersion in the dye. Real staining protocol depends on personal requests and priorities.

Preparing the sample and diagnostics

Use only appropriate instruments for collecting and preparing the samples. Process the samples with modern technology and mark them clearly. Follow the manufacturer's instructions for handling. In order to avoid mistakes, the staining procedure and diagnostics should only be conducted by authorized and qualified personnel. Use only microscope according to standards of the medical diagnostic laboratory.

Safety at work and environmental protection

Handle the product in accordance with safety at work and environmental protection guidelines. Used solutions and out of date solutions should be disposed of as special waste in accordance with national guidelines. Reagents used in this procedure could pose danger to human health. Tested tissue specimens are potentially infectious. Necessary safety measures for protecting human health should be taken in accordance with good laboratory practice. Act in accordance with signs and warnings notices printed on the product's label, as well as in BioGnost's material safety data sheet.














Storing, stability and expiry date

Keep Leishman's solution in a tightly closed original package at temperature between 15°C and 25°C. Do not keep in cold places, do not freeze and avoid exposing to direct sunlight. Date of manufacture and expiry date are printed on the product's label.

References

1. Beck, R.C. (1938): *Laboratory Manual of Hematological Technique*, Philadelphia, W.B. Saunders & Co.
2. Dacie, J. et Lewis S. (1995): *Practical haematology*, 4th ed., London, Churchill Livingstone.
3. Garcia, L. S. (2001): *Diagnostic Medical Parasitology*, 4th ed., Washington, D.C., ASM Press.
4. Giemsa, G. (1922): Das Wesen der Giemsa-Färbung, *Zentralbl f Bakt*; 89, p 99-106.
5. Kiernan, J.A. (2008): *Histological and histochemical methods: Theory and Practice*, 4th ed., Bloxham, Scion Publishing Ltd.

LE-OT-X, V10-EN8, 21 May 2019, IŠP/VR

	Refer to the supplied documentation		Storage temperature range		Number of tests in package		Product code		European Conformity
	Refer to supplied instructions		Keep away from heat and sunlight		Valid until		Lot number		Manufacturer
	For <i>in vitro</i> diagnostic use only		Keep in dry place		Caution - fragile				

 BIOGNOST Ltd.
Medjugorska 59
10040 Zagreb
CROATIA
www.biognost.com

