

SULFOSALICYLIC ACID, 20% SOLUTION

IVD *In vitro* diagnostic medical device CE

Solution for detecting proteins INSTRUCTIONS FOR USE

REF Product code: SSK-OT-100 (100 mL)

Introduction

Sulfosalicylic acid, 20% solution is most commonly used for detecting proteins in urine. The acid causes precipitation of dissolved proteins; their amount is detected by level of urine clouding.

Product description

- **SULFOSALICYLIC ACID, 20% SOLUTION** - Solution for detecting proteins in urine.
- Chemical formula: C₇H₆O₆S
- CAS: 5965-83-3

Physical properties

- Appearance: colorless to light pink fluid
- Odor: slightly pungent
- pH value (20°C): below 0.5
- Density: 1.08 g/cm³

Safety at work and environmental protection

Handle the product in accordance with safety at work and environmental protection guidelines. Used solutions and out of date solutions should be disposed of as special waste in accordance with national guidelines. Chemicals used in this procedure could pose danger to human health. Tested tissue specimens are potentially infectious. Necessary safety measures for protecting human health should be taken in accordance with good laboratory practice. Act in accordance with signs and warnings notices printed on the product's label, as well as in BioGnost's material safety data sheet.
















Storing, stability and expiry date

Keep Sulfosalicylic acid, 20% solution in a tightly closed original package at temperature between 15°C and 25°C. Do not keep in cold places, do not freeze and avoid exposing to direct sunlight. Date of manufacture and expiry date are printed on the product's label.

References

1. Orhanović, Z (2003): Kemijske štetnosti, Zagreb: Iproz,

SSK-OT-100, V2-EN2, 13 December 2018, AK/IŠP

 Refer to the supplied documentation	 Storage temperature range	 Number of tests in package	 Product code	 European Conformity	 BIOGNOST Ltd. Medjugorska 59 10040 Zagreb CROATIA www.biognost.com	
 Refer to supplied instructions	 Keep away from heat and sunlight	 Valid until	 Lot number	 Manufacturer		
 For <i>in vitro</i> diagnostic use only	 Keep in dry place	 Caution - fragile				