

BIO-DIFF KIT



IVD In vitro diagnostic medical device

Classified acc. to Regulation (EU) 2017/746 - Class A device

Three-reagent kit for rapid staining in hematology

Contains a fixative, a red, and a blue component for fast and effective polychromatic staining

INSTRUCTIONS FOR USE

BASIC UDI number							385889212HPC3010302HMCA								
EMDN code							W0103010302								
REF		Catalogue number		Volume		UDI-DI number		REF		Catalogue number		Volume		UDI-DI number	
		BD-K-100		3 x 100 mL		03858888822019				BD-K-1L		3 x 1000 mL		03858888822033	
		BD-K-500		3 x 500 mL		03858888822026				BD-K-2.5L		3 x 2500 mL		03858892120989	



Intended use and test principle

Polychromatic Romanowsky dyes are a standard in staining blood smears and bone marrow in hematology. Various sorts of Romanowsky dyes (Giemsa, May-Gruenwald, Leishman, Wright, Jenner and others) contain different ratios of methylene blue used as cation component (as well as thiazine dyes, such as azure B) and eosin Y as anion component. Cation and anion components interaction creates a well-known Romanowsky effect that cannot be achieved if each component is used individually; color purple indicates the effect's presence. Staining intensity depends on the azure B content, as well as azure B to eosin Y ratio, while a few other factors affect the result of staining: working solution and buffer solution pH value, buffer substance type, fixing method, and incubation time. BioGnost's Bio-Diff kit stains hematological samples in a short period of time and provides precise staining results, such as results of May-Gruenwald Giemsa method. Each component of the kit is stabilized and prepared according to the highest standards. Besides Bio-Diff reagents, the kit also contains buffer tablets with 6.8 and 7.2 pH values.

Product description

- BIO-DIFF KIT – Fast and efficient hematological samples staining kit.

Kit contains:	BD-K-100 (3 x 100 mL)	BD-K-500 (3 x 500 mL)	BD-K-1L (3 x 1000 mL)	BD-K-2.5L (3 x 2500 mL)
Bio-Diff 1 reagent	100 mL (BD1-OT-100)	500 mL (BD1-OT-500)	1000 mL (BD1-OT-1L)	2500 mL (BD1-OT-100)
Bio-Diff 2 reagent	100 mL (BD2-OT-100)	500 mL (BD2-OT-500)	1000 mL (BD2-OT-1L)	2500 mL (BD2-OT-100)
Bio-Diff 3 reagent	100 mL (BD3-OT-100)	500 mL (BD3-OT-500)	1000 mL (BD3-OT-1L)	2500 mL (BD3-OT-100)
Buffer tablet pH 6.8	2 pcs	5 pcs	10 pcs	15 pcs
Buffer tablet pH 7.2	2 pcs	5 pcs	10 pcs	15 pcs

Additional reagents and materials that can be used in staining

- BioGnost's immersion oils, such as BioGnost's Immersion oil, Immersion oils types C, A, FF, 37 and Tropical Grade
- VitroGnost slides and coverslips for use in histopathology and cytology

Preparation of solutions

- Buffer solution, pH 6.8 or 7.2
Dissolve 1 pH buffer tablet in 1 liter of distilled/demineralized water while stirring. Filter the solution.

Note

Pour the reagents into glass staining jars (type Coplin, Hellendahl or Schifferdecker) and return them to the original bottles after staining. Close well.

Blood smear/bone marrow sample staining procedure

1.	Let the smear dry Note: Prepare the peripheral blood smear by draining blood from a fresh blood sample	
2.	Dip the smear into Bio-Diff 1 reagent	5 x 1 second
3.	Decant the excessive reagent from the smear onto filter paper	
4.	Dip the smear into Bio-Diff 2 reagent Note: extend the incubation period if a stronger hue of red/purple is required	3 x 1 second up to 5 x 1 second
5.	Decant the excessive reagent from the smear onto filter paper	
6.	Dip the smear into Bio-Diff 3 reagent Note: decrease the incubation period if a stronger hue of red/purple is required	6 x 1 second 5 x 1 second
7.	Rinse the smear in pH 6.8 buffer solution	1 min (with agitation)
8.	Dry the smear	

Staining method for parasitology (*Leishmania*, *Toxoplasma*, *Microsporadia*) and microbiology samples (*Cryptosporidium*, *Pneumocystis carinii*)

1.	Dip the smear into Bio-Diff 1 reagent	1 min
2.	Decant the excessive reagent from the smear onto filter paper	
3.	Dip the smear into Bio-Diff 2 reagent	25 seconds
4.	Decant the excessive reagent from the smear onto filter paper	
5.	Dip the smear into Bio-Diff 3 reagent	25 seconds
6.	Rinse the smear in pH 7.2 buffer solution	1 min (with agitation)
7.	Dry the smear	

Sperm staining procedure

Preparing the sperm smear: Add 15 µL of fresh sperm sample on one side of the glass slide and create a thin and homogeneous smear. Let the smear dry (at least 10 minutes).

1.	Dip the smear into Bio-Diff 1 reagent	5 x 1 second
2.	Decant the excessive reagent from the smear onto filter paper	
3.	Dip the smear into Bio-Diff 2 reagent	5 x 1 second
4.	Decant the excessive reagent from the smear onto filter paper	
5.	Dip the smear into Bio-Diff 3 reagent	5 x 1 second
6.	Rinse the smear in pH 7.2 buffer solution	1 min (with agitation)
7.	Dry the smear	

In order to create a permanent sample, apply appropriate type of DPX medium on both stained and dried smear (BioMount DPX medium for covering/mounting cover slides). Cover the slide with VitroGnost cover glass.

Results for sperm staining

Head - homogeneous dark purple
Acrosome - light purple
Mid piece and tail - dark purple
Background - light pink

Histological sections staining procedure

a) preparation of histological sections for staining

- Fix (Formaldehyde NB 4%, Formaldehyde NB 10%) and process the tissue sample
- Embed the tissue in a paraffin block (BioWax 52/54, BioWax 56/58, BioWax Plus 56/58, BioWax Blue)
- Cut the paraffin block into 4-6 µm thin slices and mount on a VitroGnost microscope slide

b) sample staining procedure**Note: do not use Bio-Diff 1 reagent (used as fixative for non-histology samples)**

1.	Deparaffinize the section in xylene (BioClear) or xylene substitute (BioClear New)	3 exchanges, 10 min each
2.	Rehydrate in 100% alcohol (Histanol 100)	2 exchanges, 5 and 3 min
3.	Rehydrate in 95% alcohol (Histanol 95)	2 min
4.	Rehydrate in distilled/demineralized water	2 min
5.	Dip the section into Bio-Diff 2 reagent and gently stir	7 seconds
6.	Dip the section into Bio-Diff 3 reagent and gently stir	5 seconds
7.	Rinse in Buffer solution pH 7.2	1 min (with agitation)
8.	Decant the excessive reagent from the section onto filter paper	
9.	Dehydrate and differentiate in 95% alcohol (Histanol 95) while gently stirring	10 seconds
10.	Dehydrate in 100% alcohol (Histanol 100)	1 min
11.	Clear in xylene (BioClear) or xylene substitute (BioClear New)	2 exchanges, 5 min each

Immediately after clearing, apply an appropriate BioMount covering/mounting medium. If BioClear xylene was used, use one of BioGnost's xylene-based mountants (BioMount, BioMount High, BioMount M, BioMount DPX, BioMount C, or universal BioMount New). If BioClear New xylene substitute was used, the appropriate mountant is BioMount New. Cover the section with a VitroGnost cover glass.

Cytobacteriology samples staining procedure (urine, punctates, CSF)

1.	Let the cytology smear dry	
2.	Dip the smear into Bio-Diff 1 reagent	5 seconds
	Note: Incubate CSF for a longer period of time	1 min
3.	Decant the excessive reagent from the smear onto filter paper	
4.	Dip the smear into Bio-Diff 2 reagent	3 x 1 seconds (CSF 2 x 1 second)
	Note: extend the incubation period if a stronger hue of red/purple is required	do 5 x 1 second
5.	Decant the excessive reagent from the smear onto filter paper	
6.	Dip the smear into Bio-Diff 3 reagent	6 x 1 seconds (CSF 2 x 1 second)
	Note: decrease the incubation period if a stronger hue of red/purple is required	5 x 1 second
7.	Rinse the smear in pH 7.2 buffer solution	1 min (with agitation)
8.	Dry the smear	

Results for staining of cytological and histological samples

Nuclei - red to purple

Lymphocytes - plasma is colored blue

Monocytes - plasma is colored grey-blue

Neutrophil granulocytes - light purple

Eosinophil granulocytes - bright red to red-brown

Basophil granulocytes - dark purple to black

Thrombocytes - purple

Erythrocytes - reddish

Blood parasites - red (nuclei), blue (cytoplasm)

Helicobacter pylori – dark blue

Tissue cellular elements – blue to pink

Limitations

This product is intended for professional laboratory use for diagnostic purposes only. Deviations from the staining procedure described in this Instruction for use may cause differences in staining results.

Sample preparation and diagnostics

Use only appropriate instruments for collecting and preparing the samples. Process the samples using modern technology and mark them clearly. Be sure to follow the manufacturer's handling instructions. To avoid errors, staining, mounting of the slides, and diagnosis can only be carried out by qualified personnel. Use a microscope equipped according to medical diagnostic laboratory standards.

If a serious incident occurs during use of this product or as a result of its use, please report it to the manufacturer or authorized representative and competent authority.

Safety at work and environmental protection

Handle the product in accordance with occupational health and environmental protection guidelines. Used and expired solutions must be disposed of as special waste following national guidelines. Reagents used in this procedure can pose a danger to human health. The examined tissue samples are potentially infectious, and it is necessary to take the measures needed to protect human health in accordance with the guidelines of good laboratory practice. It is mandatory to read and act according to the information and warning signs printed on the product label and in the Safety Data Sheet, which is available on request.

Storage, stability, and shelf life

Upon receipt, store the product in a dry place and well-closed original packaging at a temperature of +15 °C to +25 °C. Do not freeze or expose to direct sunlight. After first opening, the product can be used until the specified expiry date, if stored properly. The expiration date is printed on the product label.

Literature

- Beck, R.C. (1938): *Laboratory Manual of Hematological Technique*, Philadelphia, W.B. Saunders & Co.
- Dacie, J. et Lewis S. (1995): *Practical haematology*, 4th ed., London, Churchill Livingstone.
- Giemsa, G. (1922): Das Wesen der Giemsa-Färbung, *Zentralb f Bakt*; 89, p 99-106.
- Kiernan, J.A. (2008): *Histological and histochemical methods: Theory and Practice*, 4th ed., Bloxham, Scion Publishing Ltd.
- May, R. et Grünwald L. (1909): *Über die Färbung von Feuchtpreparaten mit meiner Azur-Eosine methode*, Deutsche med Xschr, 35, p 1751-1752.

Warnings and precautions regarding the materials contained in the product:

	H225	Highly flammable liquid and vapour.
	H301+H311+H331	Toxic if swallowed, in contact with skin or if inhaled.
	H370	Causes damage to organs (eyes).
	P210	Keep away from heat, hot surfaces, sparks, open flames and other sources of ignition. No smoking.
	P233	Keep container tightly closed.
	P280	Wear protective gloves/protective clothing/eye protection/face protection.
	P301+P310	IF SWALLOWED: Immediately call a POISON CENTER/doctor.
	P302+P352	IF ON SKIN: Wash with plenty of soap and water.
	P304+P340	IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
	P308+P311	IF exposed or concerned: Call a POISON CENTER or doctor.

BD-K-X-IFU_EN10, 30 October 2025, IŠP

	Manufacturer		Batch code		Consult instructions for use		European conformity
	Date of manufacture		Catalogue number		Caution		Unique device identifier
	Use-by date		Temperature limit		In vitro diagnostic medical device		

 **BioGnost Ltd.**
Medjugorska 59, 10040 Zagreb, Croatia, EU, www.biognost.com

Version	Description / reason for change	Date
10.	Addition of: the QR code, table with the Warnings and precautions, revised table with the symbols, staining results of the Helicobacter pylori	30.10.2025.