

BioFix GL

FORMALIN SUBSTITUTE FOR HISTOLOGY

BioFix GL (a formalin substitute) is a next-generation reagent designed for superior tissue preservation at both structural and molecular levels. Its cutting-edge formulation ensures excellent fixation quality **while being completely non-toxic and free from any carcinogenic risks** – a safe and effective alternative to traditional fixatives like formalin. By using BioFix GL, you simultaneously enhance the protection of healthcare workers and the environment without compromising routine workflows and diagnostics.

BioFix GL

500 mL	BFGL-500
1000 mL	BFGL-1L
2500 mL	BFGL-2.5L
5000 mL	BFGL-5L

- ▶ **Effectiveness:** BioFix GL provides high-quality tissue fixation suitable for diagnostic procedures involving standard and special staining, immunohistochemistry, and molecular cytogenetic analyses
- ▶ **Versatility:** Suitable for the fixation of all tissue types
- ▶ **Compatibility:** Compatible with all processing and staining protocols
- ▶ **User-friendly:** Safe and ready to use

The shelf life of BioFix GL is 12 months when kept between +2 and +8 °C. The product can be transported at room or ambient temperature.

